



The
**Patent
Office**

GB 99/2425

4


PCT/GB99/02425



INVESTOR IN PEOPLE

09/762163

The Patent Office
Concept House
Cardiff Road
Newport
South Wales
NP10 8QQ



I, the undersigned, being an officer duly authorised in accordance with Section 74(1) and (4) of the Deregulation & Contracting Out Act 1994, to sign and issue certificates on behalf of the Comptroller-General, hereby certify that annexed hereto is a true copy of the documents as originally filed in connection with the patent application identified therein.


In accordance with the Patents (Companies Re-registration) Rules 1982, if a company named in this certificate and any accompanying documents has re-registered under the Companies Act 1980 with the same name as that with which it was registered immediately before re-registration save for the substitution as, or inclusion as, the last part of the name of the words "public limited company" or their equivalents in Welsh, references to the name of the company in this certificate and any accompanying documents shall be treated as references to the name with which it is so re-registered.

In accordance with the rules, the words "public limited company" may be replaced by p.l.c., plc, P.L.C. or PLC.

Re-registration under the Companies Act does not constitute a new legal entity but merely subjects the company to certain additional company law rules.

**PRIORITY
DOCUMENT**

SUBMITTED OR TRANSMITTED IN
COMPLIANCE WITH RULE 17.1(a) OR (b)



Signed *Andrew Gersy*

Dated 12 August 1999

This Page Blank (uspto)

The
Patent
Office

1/77

The Patent Office
Cardiff Road
Newport
Gwent NP9 1RH**Request for grant of a patent**

(See the notes on the back of this form. You can also get an explanatory leaflet from the Patent Office to help you fill in this form)

1. Your reference	80.68691		
2. Patent application number (The Patent Office will fill in this,)	9816944.4 <small>05AUG98 03000007-1 000007 001/7700 05.00 - 9816944.4</small>		
3. Full name, address and postcode of the or of each applicant (<i>underline all surnames</i>)	Dynex Technologies (Guernsey) Ltd. Rue du Commerce, St. Peter Port Guernsey, Channel Islands GY1 2HX		
Patents ADP number (<i>if you know it</i>)			
If the applicant is a corporate body, give country/state of incorporation	7488570001		
4. Title of the invention	Microplate analyser		
5. Name of your agent (<i>if you have one</i>)	Frank B. Dehn & Co.		
"Address for service" in the United Kingdom to which all correspondence should be sent (including the postcode)	179 Queen Victoria Street London EC4V 4EL		
Patents ADP number (<i>if you know it</i>)	166001		
6. If you are declaring priority from one or more earlier patent applications, give the country and the date of filing of the or of each of these earlier applications and (<i>if you know it</i>) the or each application number	Country	Priority application number (<i>if you know it</i>)	Date of filing (day / month / year)
7. If this application is divided or otherwise derived from an earlier UK application, give the number and the filing date of the earlier application	Number of earlier application	Date of filing (day / month / year)	
8. Is a statement of inventorship and of right to grant of a patent required in support of this request? (<i>Answer 'Yes' if:</i> a) any applicant named in part 3 is not an inventor, or b) there is an inventor who is not named as an applicant, or c) any named applicant is a corporate body. See note (d))	Yes		

Patents Form 1/77

9. Enter the number of sheets for any of the following items you are filing with this form. Do not count copies of the same document

Continuation sheets of this form

Description

2

Claim(s)

Abstract

Drawing(s)

1 + 1

10. If you are also filing any of the following, state how many against each item.

Priority documents

Translations of priority documents

Statement of inventorship and right to grant of a patent (*Patents Form 7/77*)

Request for preliminary examination and search (*Patents Form 9/77*)

Request for substantive examination (*Patents Form 10/77*)

Any other documents
(please specify)

11.

I/We request the grant of a patent on the basis of this application.

Signature

Date 4 August 1998

12. Name and daytime telephone number of person to contact in the United Kingdom

K.J. Tomlinson
01273 244200

Warning

After an application for a patent has been filed, the Comptroller of the Patent Office will consider whether publication or communication of the invention should be prohibited or restricted under Section 22 of the Patents Act 1977. You will be informed if it is necessary to prohibit or restrict your invention in this way. Furthermore, if you live in the United Kingdom, Section 23 of the Patents Act 1977 stops you from applying for a patent abroad without first getting written permission from the Patent Office unless an application has been filed at least 6 weeks beforehand in the United Kingdom for a patent for the same invention and either no direction prohibiting publication or communication has been given, or any such direction has been revoked.

Notes

- If you need help to fill in this form or you have any questions, please contact the Patent Office on 0645 500505.
- Write your answers in capital letters using black ink or you may type them.
- If there is not enough space for all the relevant details on any part of this form, please continue on a separate sheet of paper and write "see continuation sheet" in the relevant part(s) of the form. Any continuation sheet should be attached to this form.
- If you have answered 'Yes', *Patents Form 7/77* will need to be filed.
- Once you have filled in the form you must remember to sign and date it.
- For details of the fee and ways to pay please contact the Patent Office.

68691.345

Microplate Analyser

5 Microplate analysers are used in automated immunoassays and typically carry out automatically the sample distribution, incubation, reagent addition, washing and detection phases of immunoassays in clinical laboratories, for in-vitro diagnostic assays and other research and industrial applications.

10 Automated microplate analysers traditionally require two separate mechanisms to transport microplates between process stages (incubation, washing etc.) and dispense samples and reagents into microplates.

15 According to the invention, we provide a microplate analyser having a movable carriage carrying a pipette mechanism for dispensing fluids, and also having a mechanism allowing the carriage to pick up a microplate.
20 The latter mechanism may be any suitable clamp, engagement member or latch capable of reliably engaging, lifting and carrying a microplate. Preferably a clamp mechanism is employed.

25 Our invention allows the pipette carriage to pick up microplates, thereby saving the cost of an additional motion system.


30 An example of the invention is shown in the drawing. A pipette carriage 1 movable by suitable drives (not shown) in x, y and z directions carries a first linear actuator 2 for driving a syringe and a second linear actuator 3 for carrying a disposable dispensing tip on a conical support 4.

35 A frame 5 for carrying a microplate is provided with a mounting mechanism 6 by means of which the frame and

microplate can be picked up. In this embodiment a slot
7 is formed in the mounting mechanism and a reduced-
diameter portion 8 can enter the slot and grip the
mounting mechanism against the top of the conical
5 support 4 on retraction of the second linear actuator.

In use, the carriage 1 can be used to move around a
plate of samples as described. In addition, the
carriage can move to a supply of disposable tips, pick
10 up a tip by lowering the conical support 4 into a tip,
and dispense and aspirate fluids through the tip by
operating the syringe 2 (which is connected by a pipe
(not shown) to an axial passage through the interior of
the second actuator 3). Additionally, the second
15 actuator 3 can be positioned over a tip waste container,
and by more fully retracting the actuator the tip can be
knocked off the support 4 to fall into the waste
container.

20 If desired, only a single actuator could be used for the
three functions of syringing, plate clamping and tip
ejection, but at the expense of an additional dead
volume in the syringe pump.



PCT/SB9: 2425

Frank B. Dehn

13/8/9.

This Page Blank (uspto)